

Understanding Stability of Cys³⁴-Adducts from Intraclass Correlation Coefficient

Yeunook Bae¹, Nathan T. Montgomery¹, and William E. Funk¹


¹Department of Preventive Medicine, Feinberg School of Medicine, Northwestern University, Chicago, IL, 60611

Motivation



- Cys³⁴-adducts in the blood can be used as biomarkers for testing potential adverse health effects and environmental exposures.
- Because repeat blood sampling is not possible, blood samplings can be taken only limited numbers from patients.
- It is important to understand the stability of the Cys³⁴-adducts' levels in repeat measurements from patients during a long-term period.
- To understand the stability, intraclass correlation coefficient (ICC) value for each Cys³⁴-adduct was obtained.

Results & Discussion

Adducts with High ICC



Adducts with Low ICC




ICC Calculation & Experimental Methods

What is ICC and How to Calculate?

- Between-subject variance (σ_b^2) and within-subject variance (σ_w^2)
 - ICC close to 0: Repeated measurements from a certain subject tends to take any values across the range of distribution.
 - ICC close to 1: Repeated measurements tend to stay similar values in particular subjects.
- ICC calculation: Analysis of variance (ANOVA) by Microsoft Excel.

$$ICC = \frac{\sigma_b^2}{\sigma_b^2 + \sigma_w^2}$$



Sample Collection

- Plasma extracted from venous blood samples.
- Samples taken from 10 female volunteers: from each volunteer over a year during 2018.

Sample Pretreatment & Analysis

- Pretreat plasma with methanol to extract human serum albumin (HSA).
- Digest the HSA with trypsin.
- Analyze Cys³⁴-adducts using triple-quadrupole mass spectrometry (QQQ-MS).

- Cysteine-containing small thiol Cys³⁴-adducts (e.g., adducts of cysteine, γ -glutamylcysteine, or glutathione) or Cys³⁴-homocysteine have low ICC values.
- Adducts with low ICC have a narrow range of variation.
- Since the average among participants is similar, the samples would become good control group in terms of these adducts.
- e.g., High adduct levels from certain subjects can clearly distinguish from those of other subjects.

Conclusion

- This study used repeat sampling over a year to evaluate within- vs between-person variance of individual protein adducts using ICC.
- Several important Cys³⁴-adducts had high ICC values, suggesting single-spot measure provide robust estimates of exposure when repeat sampling is not possible.
- To understand average and population variability of certain adduct level of control group, low ICC can become a valuable tool.

Acknowledgements

We would like to thank Dr. Dale Sandler and Alexandra White at National Institute for Environmental Health Sciences for providing the blood samples. This work was supported through Grant 1 R21 ES026776-01 from the National Institute of Health.

What is Cys³⁴-Adducts?

- Cys³⁴ is one of the amino acids which are the parts of protein (peptide bond) in the human blood.
- Due to the thiol side chain, Cys³⁴ is strongly nucleophilic, which is highly reactive with the toxic electrophiles, such as environmental contaminants.
- By reacting with electrophiles, Cys³⁴ can transform into adducts, which is Cys³⁴-adducts!